

# QED 6 Plug-In



The Ansur QED 6 Plug-In is an add-on module for the Ansur Test Automation software application. The plug-in is used in conjunction with the QED 6 Defibrillator Analyzer to provide test automation capabilities using the Ansur PC software.

## Introduction

This read me file contains the latest information for Fluke Biomedical's Ansur QED 6 Plug-In.

### Serial Cable Requirements

The QED 6 requires a custom serial cable to communicate with a PC. (Part number 3370512 as of this writing.) This cable requires that a specific end of the cable be attached to the PC. Keep the orientation of the cable in mind when troubleshooting communications problems with the QED 6.

## History

This chapter describes the changes made between each new released version.

### Release 1.0.3

- Plug-in works with Laptop + Docking station.

### Release 1.0.2

- Communication problem resolved related to USB to serial converter.

### Release 1.0.1

- Fixed a problem related execution of test with different models of QED6.

### Release 1.0.0

- V1.0.0 is the first release of the QED 6 plug-in.

## Language Changes

This chapter describes the changes made to the language file in each new version. The list includes only versions where the language file was changed.

### Release 1.0.0

- Initial revision.

#### Fluke Biomedical

PO Box 9090  
Everett, WA 98206-9090 USA

Fluke Biomedical France  
30, rue Paul Claudel  
91000 EVRY, France

For more information call:  
U.S.A. (800) 648-7952 or Fax (425) 446-5629  
Europe / M-East / Africa +33 1 60 78 88 99 or Fax +33 1 60 78 68 39  
Other countries +1 (425) 446-6945 or Fax +1 (425) 446-5629  
Internet: <http://www.flukebiomedical.com>

© 2008 Fluke Biomedical. All rights reserved.