

# ESA601 Plug-In Read Me

**FLUKE.**

Biomedical

Ansur ESA601 Plug-In is a windows application facilitating automatic safety testing via the ESA601 Electrical Safety Analyzer.

## Introduction

This read me file contains last minute information for Fluke Biomedical's Ansur ESA601 Plug-In.

## History

This chapter describes the changes made between each new released version.

### Release 1.1.0

- Upgraded the test guide to support the new test guide toolbar introduced with Ansur V2.4.
- The plug-in now supports COM ports with a port number greater than 16.

### Release 1.0.3

- Upgraded to support changes in firmware version 1.05 related to MAP test.

### Release 1.0.2

- Fixed a problem in communication where patient lead 2 and 3 was interchanged.

### Release 1.0.1

- In Test Lead Calibration, the Calibrate button is now hidden once a calibration has been successfully completed.
- Added sample templates for the different standards provided with the plug-in.
- User manual was improved.

### Release 1.0.0

- V1.0.0 is the first release of the ESA601 plug-in.

## Language Changes

This chapter describes the changes made to the language file in each new version. The list includes only versions where the language file was changed.

### Release 1.1.0

- Altered: 20307

#### Fluke Corporation

PO Box 9090  
Everett, WA USA 98206

Fluke Eurpoe B.V.  
PO Box 1186, 5602 BD  
Eindhoven, The Netherlands

For more information call:  
U.S.A. (800) 443-5853 or Fax (425) 356-5116  
Europe / M-East (31 40) 2 678 200 or Fax (31 40) 2 678 222  
Canada (905) 890-7600 or Fax (905) 890-6866  
Other countries (425) 356-5500 or Fax (425) 356-5116  
Internet: <http://www.fluke.com>